

Алматы (7273)495-231  
Ангарск (3955)60-70-56  
Архангельск (8182)63-90-72  
Астрахань (8512)99-46-04  
Барнаул (3852)73-04-60  
Белгород (4722)40-23-64  
Благовещенск (4162)22-76-07  
Брянск (4832)59-03-52  
Владивосток (423)249-28-31  
Владикавказ (8672)28-90-48  
Владимир (4922) 49-43-18  
Волгоград (844)278-03-48  
Вологда (8172)26-41-59  
Воронеж (473)204-51-73  
Екатеринбург (343)384-55-89

Ижевск (3412)26-03-58  
Иваново (4932)77-34-06  
Иркутск (395)279-98-46  
Казань (843)206-01-48  
Калининград (4012)72-03-81  
Калуга (4842)92-23-67  
Кемерово (3842)65-04-62  
Киров (8332)68-02-04  
Коломна (4966)23-41-49  
Кострома (4942)77-07-48  
Краснодар (861)203-40-90  
Красноярск (391)204-63-61  
Курск (4712)77-13-04  
Курган (3522)50-90-47  
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13  
Москва (495)268-04-70  
Мурманск (8152)59-64-93  
Набережные Челны (8552)20-53-41  
Нижний Новгород (831)429-08-12  
Новокузнецк (3843)20-46-81  
Ноябрьск (3496)41-32-12  
Новосибирск (383)227-86-73  
Ноябрьск (3496)41-32-12  
Омск (3812)21-46-40  
Орел (4862)44-53-42  
Оренбург (3532)37-68-04  
Пенза (8412)22-31-16  
Петрозаводск (8142)55-98-37  
Псков (8112)59-10-37

Пермь (342)205-81-47  
Ростов-на-Дону (863)308-18-15  
Рязань (4912)46-61-64  
Самара (846)206-03-16  
Саранск (8342)22-96-24  
Санкт-Петербург (812)309-46-40  
Саратов (845)249-38-78  
Севастополь (8692)22-31-93  
Симферополь (3652)67-13-56  
Смоленск (4812)29-41-54  
Сочи (862)225-72-31  
Ставрополь (8652)20-65-13  
Сыктывкар (8212)25-95-17  
Сургут (3462)77-98-35  
Тамбов (4752)50-40-97

Тверь (4822)63-31-35  
Тольяти (8482)63-91-07  
Томск (3822)98-41-53  
Тула (4872)33-79-87  
Тюмень (3452)66-21-18  
Улан-Удэ (3012)59-97-51  
Ульяновск (8422)24-23-59  
Уфа (347)229-48-12  
Хабаровск (4212)92-98-04  
Чебоксары (8352)28-53-07  
Челябинск (351)202-03-61  
Череповец (8202)49-02-64  
Чита (3022)38-34-83  
Якутск (4112)23-90-97  
Ярославль (4852)69-52-93

Киргизия (996)312-96-26-47

Россия (495)268-04-70

Казахстан (772)734-952-31

<https://supelco.nt-rt.ru> || [suz@nt-rt.ru](mailto:suz@nt-rt.ru)

## Капиллярные неполярные колонки



**Non-polar GC columns** are best suited for the analysis of non-polar compounds (e.g. alkanes). Equity<sup>®</sup>-1, Omegawax<sup>®</sup>, Petrocol<sup>®</sup>, SPB<sup>®</sup>-1, SPB<sup>®</sup>-Octyl, SLB<sup>®</sup>-5ms, SUPELCO WAX<sup>®</sup> and VOCOL<sup>®</sup> are highly recommended products for the separation of non-polar compounds. Polar GC columns are most efficient at separating polar analytes. SPB<sup>®</sup>-624, OVI-G43, VOCOL<sup>®</sup>, SPB<sup>®</sup>-20, Equity<sup>®</sup>-1701, SPB<sup>®</sup>-35, SPB<sup>®</sup>-50, SPB<sup>®</sup>-225, PAG, Omegawax<sup>®</sup>, SPB<sup>®</sup>-1000, SLB<sup>®</sup>-IL60i, Nukol<sup>™</sup>, and SUPELCO WAX<sup>®</sup> 10 are some of the polar columns that provide the best separation of polar compounds.

Highly polar columns include those that contain stationary phases with a higher percentage of cyanopropyl functional groups. These phases exhibit a variety of interactions amenable to the separation of polarizable compounds. Polarizable compounds are those containing C and H, but with one or more C=C or C≡C bonds. Highly polar capillary columns include SP<sup>®</sup>-2330, SP<sup>®</sup>-2331, SP<sup>®</sup>-2380, SP<sup>®</sup>-2560, SP<sup>®</sup>-2340, TCEP, SLB<sup>®</sup>-IL1111i.

## Non-Polar Fused Silica Tubing

Compare Product No.	Description
<input type="checkbox"/> 25758	L × I.D. 15 m × 0.53 mm
<input type="checkbox"/> 25743	L × I.D. 5.0 m × 0.32 mm

Compare Product No.	Description
<input type="checkbox"/> 25307	L × I.D. 1 m × 0.53 mm
<input type="checkbox"/> 25744	L × I.D. 5.0 m × 0.53 mm
<input type="checkbox"/> 25742	L × I.D. 5.0 m × 0.25 mm



## Polar Fused Silica Tubing

Compare Product No.	Description
<input type="checkbox"/> 25777	L × I.D. 30 m × 0.25 mm
<input type="checkbox"/> 25765	L × I.D. 15 m × 0.32 mm
<input type="checkbox"/> 25710	L × I.D. 1 m × 0.10 mm
<input type="checkbox"/> 25714	L × I.D. 1 m × 0.53 mm
<input type="checkbox"/> 25753	L × I.D. 5.0 m × 0.53 mm

## Intermediate Polar Fused Silica Tubing

Compare Product No.	Description
<input type="checkbox"/> 25729	L × I.D. 3.0 m × 0.53 mm
<input type="checkbox"/> 25747	L × I.D. 5.0 m × 0.25 mm
<input type="checkbox"/> 25760-U	L × I.D. 15 m × 0.25 mm
<input type="checkbox"/> 25339	L × I.D. 5.0 m × 0.53 mm
<input type="checkbox"/> 25708	L × I.D. 1 m × 0.32 mm



## F.A.M.E. Column Evaluation Mix

200-838-9

Compare Product No.	Description
<input type="checkbox"/> 47088-U	certified reference material, 1000 µg/mL each component in dichloromethane



## Omegawax™ Column Test Mix

203-777-6

Compare Product No.	Description
<input type="checkbox"/> 48476	50 µg/mL in hexane, analytical standard



## SP<sup>®</sup>-2560 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 23362-U	L × I.D. 100 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 24056	L × I.D. 100 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 23348-U	L × I.D. 75 m × 0.18 mm, d <sub>f</sub> 0.14 μm
<input type="checkbox"/> 29688-U	L × I.D. 200 m × 0.25 mm, d <sub>f</sub> 0.20 μm



## Capillary Column Nut (fits Agilent injectors and non-MS detectors), hexagonal wrenchtight version

Compare Product No.	Description
<input type="checkbox"/> 24833-U	pkg of 2 ea



## Supelco Ferrule Nut Adapter (fits Agilent injectors and non-MS detectors), knurled fingertight version

Compare Product No.	Description
<input type="checkbox"/> 22509	pkg of 2 ea



## Supelco Ferrule Nut Adapter (fits Agilent injectors and non-MS detectors), hexagonal wrenchtight version

Compare Product No.	Description
<input type="checkbox"/> 22470-U	pkg of 2 ea



## SPB<sup>®</sup>-5 Intuvo Capillary GC Column

Synonym(s): SPB<sup>®</sup>-5 Intuvo Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24048INT	L × I.D. 30 m × 0.32 mm, d <sub>f</sub> 0.25 μm



## Equity<sup>®</sup>-5 Intuvo Capillary GC Column

Synonym(s): Equity<sup>®</sup>-5 Intuvo Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 28089UINT	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm, for use with Agilent Intuvo GC



## SPB<sup>®</sup>-1 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24229-U	L × I.D. 12 m × 0.20 mm, d <sub>f</sub> 0.33 μm
<input type="checkbox"/> 23302-U	L × I.D. 60 m × 0.75 mm, d <sub>f</sub> 1.00 μm
<input type="checkbox"/> 24134-U	L × I.D. 10 m × 0.20 mm, d <sub>f</sub> 1.20 μm
<input type="checkbox"/> 24162	L × I.D. 15 m × 0.20 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 24163	L × I.D. 30 m × 0.20 mm, d <sub>f</sub> 0.20 μm



## SPB<sup>®</sup>-5 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 25420-U	L × I.D. 30 m × 0.53 mm, d <sub>f</sub> 1.00 μm
<input type="checkbox"/> 24166	L × I.D. 30 m × 0.20 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 24032	L × I.D. 15 m × 0.25 mm, d <sub>f</sub> 0.25 μm
<input type="checkbox"/> 24034	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm
<input type="checkbox"/> 24036	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## Petrocol<sup>®</sup> DH Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24160-U	L × I.D. 100 m × 0.25 mm, d <sub>f</sub> 0.50 μm



## SAC™ -5 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24156	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## SP®-2330 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24020-U	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 24019	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 24018	L × I.D. 15 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 24073	L × I.D. 30 m × 0.32 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 24074	L × I.D. 60 m × 0.32 mm, d <sub>f</sub> 0.20 μm



## SP®-2340 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24022	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 24023	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 24076	L × I.D. 60 m × 0.32 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 24075	L × I.D. 30 m × 0.32 mm, d <sub>f</sub> 0.20 μm



## PAG Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24223	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## SUPELCOWAX™ 10 Intuvo Capillary GC Column

Synonym(s): SUPELCOWAX™ 10 Intuvo Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24080UINT	column L × I.D. 30 m × 0.32 mm, d <sub>f</sub> 0.25 μm, for use with Agilent Intuvo GC
<input type="checkbox"/> 24084INT	column L × I.D. 30 m × 0.32 mm, d <sub>f</sub> 0.50 μm, for use with Advanced Pharmaceutical Sciences Intuvo GC



## SP<sup>®</sup>-2380 Intuvo Capillary GC Column

Synonym(s): SP<sup>®</sup>-2380 Intuvo Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24111INT	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.20 μm, for use with Agilent Intuvo GC
<input type="checkbox"/> 24110UINT	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm, for use with Agilent Intuvo GC



## SLB<sup>®</sup>-5ms Intuvo Capillary GC Column

Synonym(s): SLB<sup>®</sup>-5ms Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 28472UINT	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.25 μm, for use with Agilent Intuvo GC
<input type="checkbox"/> 28564UINT	L × I.D. 20 m × 0.18 mm, d <sub>f</sub> 0.18 μm, for use with Agilent Intuvo GC
<input type="checkbox"/> 28482UINT	L × I.D. 30 m × 0.32 mm, d <sub>f</sub> 0.25 μm, for use with Agilent Intuvo GC
<input type="checkbox"/> 28473UINT	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.50 μm, for use with Agilent Intuvo GC



## SPB<sup>®</sup>-1 Intuvo Capillary GC Column

Synonym(s): SPB<sup>®</sup>-1 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24028INT	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm, for use with Agilent Intuvo GC



[All Photos](#)(3)

## Equity<sup>®</sup>-1 Intuvo Capillary GC Column

Synonym(s): Equity<sup>®</sup>-1 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 28046UINT	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm, for use with Agilent Intuvo GC



## SLB<sup>®</sup>-1ms Intuvo Capillary GC Column

Synonym(s): SLB<sup>®</sup>-1ms Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 29731UINT	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm, for use with Agilent Intuvo GC



## SP<sup>®</sup>-2340 Intuvo Capillary GC Column

Synonym(s): SP<sup>®</sup>-2340 Capillary GC Column

Compare Product No.	Description
---------------------	-------------

<input type="checkbox"/> 24023INT	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.20 μm, for use with Agilent Intuvo GC
-----------------------------------	---



## Astec<sup>®</sup> ChiralDEX<sup>™</sup> G-TA Intuvo Capillary GC Column

Synonym(s): Astec<sup>®</sup> CHIRALDEX<sup>™</sup> G-TA Capillary GC Column

Compare Product No.	Description
---------------------	-------------

<input type="checkbox"/> 73033ASTINT	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.12 μm, for use with Agilent Intuvo GC
--------------------------------------	---



## SLB<sup>®</sup>-IL60i Intuvo Capillary GC Column

Synonym(s): SLB<sup>®</sup>-IL60i Capillary GC Column

Compare Product No.	Description
---------------------	-------------

<input type="checkbox"/> 29833UINT	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.20 μm, for use with Agilent Intuvo GC
------------------------------------	---

<input type="checkbox"/> 29832UINT	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm, for use with Agilent Intuvo GC
------------------------------------	---



## SPB<sup>®</sup>-Octyl Intuvo Capillary GC Column

Synonym(s): SPB<sup>®</sup>-Octyl Capillary GC Column

Compare Product No.	Description
---------------------	-------------

<input type="checkbox"/> 24218UINT	column L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm, for use with Agilent Intuvo GC
------------------------------------	--



## SLB<sup>®</sup>-ILPAH Intuvo Capillary GC Column

Synonym(s): SLB<sup>®</sup>-ILPAH Capillary GC Column

Compare Product No.	Description
---------------------	-------------

<input type="checkbox"/> 29799UINT	L × I.D. 20 m × 0.18 mm, d <sub>f</sub> 0.05 μm, for use with Agilent Intuvo GC
------------------------------------	---



## Astec<sup>®</sup> CHIRALDEX<sup>™</sup> B-DM Intuvo Capillary GC Column

Synonym(s): Astec<sup>®</sup> CHIRALDEX<sup>™</sup> B-DM Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 77023ASTINT	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.12 μm, for use with Agilent (Intuvo GC)

## beta-DEX<sup>™</sup> 120 Intuvo Capillary GC Column

Synonym(s): β-DEX<sup>™</sup> 120

Compare Product No.	Description
<input type="checkbox"/> 24304INT	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm, for use with Agilent Intuvo GC



## Omegawax<sup>™</sup> Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 23399-U	L × I.D. 15 m × 0.10 mm, d <sub>f</sub> 0.10 μm

## SLB<sup>®</sup>-IL111i Intuvo Capillary GC Column

Synonym(s): SLB<sup>®</sup>-IL111i Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 29884UINT	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.20 μm, for use with Agilent Intuvo GC



## MET-SimDis Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 29894-U	L × I.D. 5 m × 0.53 mm, d <sub>f</sub> 0.88 μm
<input type="checkbox"/> 29893-U	L × I.D. 5 m × 0.53 mm, d <sub>f</sub> 0.17 μm
<input type="checkbox"/> 29895-U	L × I.D. 10 m × 0.53 mm, d <sub>f</sub> 0.88 μm
<input type="checkbox"/> 29896-U	L × I.D. 5 m × 0.53 mm, d <sub>f</sub> 2.65 μm
<input type="checkbox"/> 29898-U	L × I.D. 10 m × 0.53 mm, d <sub>f</sub> 2.65 μm



## SLB<sup>®</sup>-1ms Capillary GC Column



Compare Product No.	Description
<input type="checkbox"/> 29737-U	L × I.D. 60 m × 0.32 mm, d <sub>f</sub> 0.25 μm
<input type="checkbox"/> 29734-U	L × I.D. 30 m × 0.32 mm, d <sub>f</sub> 0.10 μm
<input type="checkbox"/> 29735-U	L × I.D. 15 m × 0.32 mm, d <sub>f</sub> 0.25 μm



## SPB<sup>®</sup>-20 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24196-U	size × I.D. 30 m × 0.25 mm, d <sub>f</sub> 1.00 μm
<input type="checkbox"/> 24086	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## SLB<sup>®</sup>-IL59 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 28880-U	L × I.D. 15 m × 0.10 mm, d <sub>f</sub> 0.08 μm
<input type="checkbox"/> 28891-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm



## SLB<sup>®</sup>-IL76i Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 29869-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm



## SLB<sup>®</sup>-IL60i Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 29832-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 29829-U	L × I.D. 20 m × 0.18 mm, d <sub>f</sub> 0.14 μm
<input type="checkbox"/> 29833-U	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 29836-U	L × I.D. 30 m × 0.32 mm, d <sub>f</sub> 0.26 μm



## SLB<sup>®</sup>-IL76 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 28913-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm



## β-DEX™ 120

231-493-2

Compare Product No.	Description
<input type="checkbox"/> 24305-U	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.25 μm
<input type="checkbox"/> 24304	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## Watercol™ 1910 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 29711-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 29714-U	L × I.D. 30 m × 0.32 mm, d <sub>f</sub> 0.26 μm



## SLB®-IL61 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 29486-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm



## SLB®-ILPAH Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 29799-U	L × I.D. 20 m × 0.18 mm, d <sub>f</sub> 0.05 μm



## SLB®-IL111i Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 29883-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 29884-U	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.20 μm



## SPB<sup>®</sup>-1 SULFUR Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24158	L × I.D. 30 m × 0.32 mm, d <sub>f</sub> 4.00 μm



## TCEP Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24161	L × I.D. 60 m × 0.32 mm, d <sub>f</sub> 0.51 μm
<input type="checkbox"/> 24153	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.44 μm



## MET-Biodiesel Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 28804-U	L × I.D. 10 m × 0.32 mm , with 2 m x 0.53 mm I.D. guard, d <sub>f</sub> 0.10 μm
<input type="checkbox"/> 28808-U	L × I.D. 15 m × 0.32 mm , with 2 m x 0.53 mm I.D. guard, d <sub>f</sub> 0.10 μm



## SLB<sup>®</sup>-IL82 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 29479-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm



## SP<sup>®</sup>-2250 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24010	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm



## SLB<sup>®</sup>-IL60 Capillary GC Column

Compare Product No.	Description
---------------------	-------------

Compare Product No.	Description
<input type="checkbox"/> 29503-U	L × I.D. 15 m × 0.10 mm, d <sub>f</sub> 0.08 μm
<input type="checkbox"/> 29504-U	L × I.D. 20 m × 0.18 mm, d <sub>f</sub> 0.14 μm
<input type="checkbox"/> 29505-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 29506-U	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 29508-U	L × I.D. 30 m × 0.32 mm, d <sub>f</sub> 0.26 μm



## SPB<sup>®</sup>-Octyl Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24233-U	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 1.00 μm
<input type="checkbox"/> 24218-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## Petrocol<sup>®</sup> EX2887 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 25337	L × I.D. 5.0 m × 0.53 mm, d <sub>f</sub> 0.10 μm



## Petrocol<sup>®</sup> 2887 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 25323	L × I.D. 5.0 m × 0.53 mm, d <sub>f</sub> 0.50 μm



## SLB<sup>®</sup>-IL100 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 28884-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 28886-U	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 28887-U	L × I.D. 30 m × 0.32 mm, d <sub>f</sub> 0.26 μm



## Astec<sup>®</sup> CHIRALDEX<sup>™</sup> GC Column Screening Kit

Compare Product No.	Description
<input type="checkbox"/> 71030AST	L × I.D. 30 m × 0.25 mm, pkg of 1 kit



## SLB<sup>®</sup>-ILD3606 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 29691-U	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.20 μm



## Petrocol<sup>®</sup> DH Octyl Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24282	L × I.D. 100 m × 0.25 mm, d <sub>f</sub> 0.50 μm



## Astec<sup>®</sup> CHIRALDEX<sup>™</sup> B-DM Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 77023AST	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.12 μm
<input type="checkbox"/> 77025AST	L × I.D. 50 m × 0.25 mm, d <sub>f</sub> 0.12 μm



## Petrocol<sup>®</sup> DH 150 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24155	L × I.D. 150 m × 0.25 mm, d <sub>f</sub> 1.00 μm

## SLB<sup>®</sup>-IL111 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 28927-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 28928-U	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 29647-U	L × I.D. 100 m × 0.25 mm, d <sub>f</sub> 0.20 μm
<input type="checkbox"/> 29689-U	L × I.D. 200 m × 0.25 mm, d <sub>f</sub> 0.20 μm



## Carbowax<sup>®</sup> Amine Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 25353	L × I.D. 30 m × 0.53 mm, d <sub>f</sub> 1.00 μm

## Petrocol<sup>®</sup> DH 50.2 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24133-U	L × I.D. 50 m × 0.20 mm, d <sub>f</sub> 0.50 μm



## α-DEX<sup>™</sup> 120

262-373-8

Compare Product No.	Description
<input type="checkbox"/> 24310	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## Astec<sup>®</sup> CHIRALDEX<sup>™</sup> B-TA Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 73023AST	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.12 μm



## Astec<sup>®</sup> CHIRALDEX<sup>™</sup> G-TA Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 73033AST	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.12 μm
<input type="checkbox"/> 73032AST	L × I.D. 20 m × 0.25 mm, d <sub>f</sub> 0.12 μm
<input type="checkbox"/> 73035AST	L × I.D. 50 m × 0.25 mm, d <sub>f</sub> 0.12 μm



## γ-DEX<sup>™</sup> 120

Compare Product No.	Description
<input type="checkbox"/> 24307	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## Astec<sup>®</sup> CHIRALDEX<sup>™</sup> A-TA Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 73003AST	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.12 μm



## β-DEX<sup>™</sup> 110

262-373-8

Compare Product No.	Description
<input type="checkbox"/> 24302	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.25 μm
<input type="checkbox"/> 24301	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## PTA-5 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24330	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 1.00 μm
<input type="checkbox"/> 24277	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.50 μm
<input type="checkbox"/> 24333	L × I.D. 30 m × 0.32 mm, d <sub>f</sub> 1.50 μm
<input type="checkbox"/> 25439	L × I.D. 30 m × 0.53 mm, d <sub>f</sub> 3.00 μm



## Astec<sup>®</sup> CHIRALDEX<sup>™</sup> B-PM Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 76023AST	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.12 μm



## Astec<sup>®</sup> CHIRALDEX<sup>™</sup> G-BP Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 75033AST	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.12 μm



## Astec<sup>®</sup> CHIRALDEX<sup>™</sup> G-PN Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 74033AST	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.12 μm



## β-DEX<sup>™</sup> 225

Compare Product No.	Description
<input type="checkbox"/> 24348	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## α-DEX<sup>™</sup> 225

262-373-8

Compare Product No.	Description
<input type="checkbox"/> 24311	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## Astec<sup>®</sup> CHIRALDEX<sup>™</sup> A-DA Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 72003AST	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.12 μm



## Astec<sup>®</sup> CHIRALDEX<sup>™</sup> B-DA Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 72023AST	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.12 μm



## Astec<sup>®</sup> CHIRALDEX<sup>™</sup> B-PH Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 71023AST	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.12 μm





## $\beta$ -DEX™ 325

Compare Product No.	Description
<input type="checkbox"/> 24308	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## $\gamma$ -DEX™ 325

Compare Product No.	Description
<input type="checkbox"/> 24306	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## SUPELCOWAX™ 10 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 24079	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## Equity®-5 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 28279-U	L × I.D. 30 m × 0.53 mm, d <sub>f</sub> 5.00 μm
<input type="checkbox"/> 28089-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm
<input type="checkbox"/> 28094-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 1.00 μm
<input type="checkbox"/> 28097-U	L × I.D. 30 m × 0.32 mm, d <sub>f</sub> 0.25 μm
<input type="checkbox"/> 28259-U	L × I.D. 30 m × 0.53 mm, d <sub>f</sub> 0.50 μm



## SLB®-5ms Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 28466-U	L × I.D. 15 m × 0.10 mm, d <sub>f</sub> 0.10 μm
<input type="checkbox"/> 28465-U	L × I.D. 10 m × 0.10 mm, d <sub>f</sub> 0.10 μm
<input type="checkbox"/> 28564-U	L × I.D. 20 m × 0.18 mm, d <sub>f</sub> 0.18 μm
<input type="checkbox"/> 28467-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.10 μm
<input type="checkbox"/> 28469-U	L × I.D. 15 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## Astec<sup>®</sup> CHIRALDEX<sup>™</sup> B-DP Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 78023AST	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.12 μm



## Astec<sup>®</sup> CHIRALDEX<sup>™</sup> G-DP Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 78033AST	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.12 μm



## γ-DEX<sup>™</sup> 225

Compare Product No.	Description
<input type="checkbox"/> 24312	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm



## Polar Column Test Mix

200-838-9

Compare Product No.	Description
<input type="checkbox"/> 47302	certified reference material, 500 μg/mL each component in dichloromethane



## Intermediate Polar Column Test Mix

200-838-9


Compare Product No.	Description
<input type="checkbox"/> 47301	certified reference material, 500 μg/mL in dichloromethane



## Equity<sup>®</sup>-1 Capillary GC Column

Compare Product No.	Description
<input type="checkbox"/> 28063-U	L × I.D. 60 m × 0.32 mm, d <sub>f</sub> 5.00 μm

Compare Product No.	Description
<input type="checkbox"/> 28041-U	L × I.D. 12 m × 0.20 mm, d <sub>f</sub> 0.33 μm
<input type="checkbox"/> 28046-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 0.25 μm
<input type="checkbox"/> 28047-U	L × I.D. 60 m × 0.25 mm, d <sub>f</sub> 0.25 μm
<input type="checkbox"/> 28049-U	L × I.D. 30 m × 0.25 mm, d <sub>f</sub> 1.00 μm



## tert-Butyl Methyl Ether

**Synonym(s):** *tert*-Butyl methyl ether, MTBE, Methyl *tert*-butyl ether

**Linear Formula:** (CH<sub>3</sub>)<sub>3</sub>COCH<sub>3</sub>

1634-04-4

88.15

216-653-1

1730942

Compare Product No.	Description
<input type="checkbox"/> PHR1560	Pharmaceutical Secondary Standard; Certified Reference Material

## Benzene

**Empirical Formula (Hill Notation):** C<sub>6</sub>H<sub>6</sub>


71-43-2

78.11

200-753-7

969212

Compare Product No.	Description
<input type="checkbox"/> PHR1310	Pharmaceutical Secondary Standard; Certified Reference Material



## L-Isoleucine

**Synonym(s):** (2S,3S)-2-Amino-3-methylpentanoic acid

**Linear Formula:** C<sub>2</sub>H<sub>5</sub>CH(CH<sub>3</sub>)CH(NH<sub>2</sub>)CO<sub>2</sub>H

73-32-5

131.17

200-798-2

1721792

Compare Product No.	Description
<input type="checkbox"/> PHR1099	Pharmaceutical Secondary Standard; Certified Reference Material

## Cyclohexane


**Synonym(s):** Cyclohexane

**Empirical Formula (Hill Notation):** C<sub>6</sub>H<sub>12</sub>

110-82-7

84.16  
203-806-2  
1900225

Compare Product No.	Description
<input type="checkbox"/> PHR1687	Pharmaceutical Secondary Standard; Certified Reference Material




## p-Xylene

Synonym(s): 1,4-Dimethylbenzene

Linear Formula: C<sub>6</sub>H<sub>4</sub>(CH<sub>3</sub>)<sub>2</sub>

106-42-3  
106.17  
203-396-5  
1901563

Compare Product No.	Description
<input type="checkbox"/> PHR1314	Pharmaceutical Secondary Standard; Certified Reference Material



## o-Xylene

Synonym(s): 1,2-Dimethylbenzene

Linear Formula: C<sub>6</sub>H<sub>4</sub>(CH<sub>3</sub>)<sub>2</sub>

95-47-6  
106.17  
202-422-2  
1815558

Compare Product No.	Description
<input type="checkbox"/> PHR1312	Pharmaceutical Secondary Standard; Certified Reference Material

Алматы (7273)495-231  
Ангарск (3955)60-70-56  
Архангельск (8182)63-90-72  
Астрахань (8512)99-46-04  
Барнаул (3852)73-04-60  
Белгород (4722)40-23-64  
Благовещенск (4162)22-76-07  
Брянск (4832)59-03-52  
Владивосток (423)249-28-31  
Владикавказ (8672)28-90-48  
Владимир (4922) 49-43-18  
Волгоград (844)278-03-48  
Вологда (8172)26-41-59  
Воронеж (473)204-51-73  
Екатеринбург (343)384-55-89

Ижевск (3412)26-03-58  
Иваново (4932)77-34-06  
Иркутск (395)279-98-46  
Казань (843)206-01-48  
Калининград (4012)72-03-81  
Калуга (4842)92-23-67  
Кемерово (3842)65-04-62  
Киров (8332)68-02-04  
Коломна (4966)23-41-49  
Кострома (4942)77-07-48  
Краснодар (861)203-40-90  
Красноярск (391)204-63-61  
Курск (4712)77-13-04  
Курган (3522)50-90-47  
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13  
Москва (495)268-04-70  
Мурманск (8152)59-64-93  
Набережные Челны (8552)20-53-41  
Новокузнецк (3843)20-46-81  
Ноябрьск (3496)41-32-12  
Новосибирск (383)227-86-73  
Ноябрьск (3496)41-32-12  
Омск (3812)21-46-40  
Орел (4862)44-53-42  
Оренбург (3532)37-68-04  
Пенза (8412)22-31-16  
Петрозаводск (8142)55-98-37  
Псков (8112)59-10-37

Пермь (342)205-81-47  
Ростов-на-Дону (863)308-18-15  
Рязань (4912)46-61-64  
Самара (846)206-03-16  
Саранск (8342)22-96-24  
Санкт-Петербург (812)309-46-40  
Саратов (845)249-38-78  
Севастополь (8692)22-31-93  
Симферополь (3652)67-13-56  
Смоленск (4812)29-41-54  
Сочи (862)225-72-31  
Ставрополь (8652)20-65-13  
Сыктывкар (8212)25-95-17  
Сургут (3462)77-98-35  
Тамбов (4752)50-40-97

Тверь (4822)63-31-35  
Тольяти (8482)63-91-07  
Томск (3822)98-41-53  
Тула (4872)33-79-87  
Тюмень (3452)66-21-18  
Улан-Удэ (3012)59-97-51  
Ульяновск (8422)24-23-59  
Уфа (347)229-48-12  
Хабаровск (4212)92-98-04  
Чебоксары (8352)28-53-07  
Челябинск (351)202-03-61  
Череповец (8202)49-02-64  
Чита (3022)38-34-83  
Якутск (4112)23-90-97  
Ярославль (4852)69-52-93

Киргизия (996)312-96-26-47

Россия (495)268-04-70

Казахстан (772)734-952-31